

Equipment required

- 5 x sterile Z serum clot tubes (gas sterilised)
- 4 x red top plain vacutainer (gas sterilised)
- 0.2-micron filters (2-4 required)
- LL Neg connector
- 3-inch blunt cannula
- 20ml syringe
- Transfer needle (any guage)
- Bench top incubator (37 °C)



Figure 1
Sterile Z serum tubes

How to: Use Z serum tubes to process ACS (autologous conditioned serum)

Introduction

In a recent study by the Norwegian School of Veterinary Science (Fjordbakk *et al.* 2014) it was shown that serum from blood incubated in plastic vacutainer tubes containing Z Serum Clot Activator contained similar levels of the cytokines and growth factors investigated as serum incubated in two specialised ACS containers. These results corroborate those of a previous study (Hraha *et al.* 2011) and indicate that cytokine production may be a consequence of whole blood incubation rather than an effect of specialised ACS containers. Z Serum clot activated tubes were capable of producing the same output as specialised ACS containers if not more.

Technique

- Using the Z serum tubes (Figure 1) collect 7.5ml of blood in each of the four tubes.
- Write the name of the patient and time of blood collection on each of the tubes and place in the incubator as soon as possible.
- Keep the tubes in an upright position while transporting and in the incubator (Figure 2).
- Incubation starts once the clot is formed. ACS preparation should be performed no later than 8hrs after the time of blood collection.
- After 8hrs remove the 5 tubes from the incubator keeping them in the tube rack (Figure 3). Place the rack on a clean, sterile surface.



Figure 2
Z serum tubes in the incubator




Figure 3
Z serum tubes in a tube rack

“Preparation of ACS must be performed no later than 8 hours after the time of blood collection”

Warning
This process involves the preparation of a substance intended for administration to a patient. Sterility and aseptic technique is paramount. Nupsala cannot be held responsible for operator error.

Processing ACS Using Z Serum Tubes

- Using the Harvest SmartPREP 2 centrifuge, use two counter weights of equal size. Place a serum tube into each of the center tube retaining ports (Figure 4). Close the lid and spin on the PRP cycle. If not using the Harvest centrifuge the g-rating will need to be calculated to give RPM. (2500g for 10min). Call the Nupsala office for help with g/RPM calculations (+44 (0)1865 922227). Once spun, remove the serum tube and place back into the rack (Figure 5).
- Using the 3-inch cannula attached to a 20ml syringe, penetrate the red rubber bung at the top of the tube, keeping the tube in an upright position (Figure 6). Extract all available serum (c2.5ml per tube – 4 tubes giving a total volume of c10ml).



Figure 4




Figure 5



Figure 6

- Place a 0.2-micron filter on the 20 ml syringe containing the serum. Place the needle on the end of the filter (Figure 7). Penetrate the top of a red plain vacutainer and transfer 2.5ml (or equal quantities /4) of the plasma, having passed through the filter, into the vacutainer (Figure 8). Do this for all four vaccutainers. Label each vaccutainer and place in a clear zip-locked bag with the patient's ID clearly labelled on the bag. Place the bag containing the tubes in the freezer.
- ACS should be used 10 days apart. Remove a vaccutainer containing serum from the freezer and thaw thoroughly (30 mins at room temp). Using a 3-inch cannula and a 5ml syringe remove the thawed serum from the red vaccutainer (Figure 9).
- Remove the cannula and place a 0.2-micron filter onto the syringe. Place a L-L connector on the remaining filter port and connect a 5ml syringe to the L-L connector. Transfer the serum from one syringe to the other via the filter (Figure 10). Remove the syringe now containing the serum and place a stopper or needle onto the syringe. The serum is now ready for administering to the patient.



Figure 7




Figure 8



Figure 9



Figure 10

References

- C. T. FJORDBAKK, G. M. JOHANSEN†, A. C. LØVÅS, K. L. OPPEGÅRD and A. K. STORSET (2014) Surgical stress influences cytokine content in autologous conditioned serum Equine Vet J.
- Hraha, T.H., Doremus, K.M., McIlwraith, C.W. and Frisbie, D.D. (2011) Autologous conditioned serum: the comparative cytokine profiles of two commercial methods (IRAP and IRAP II) using equine blood. Equine Vet. J. 43, 516-521.

Nupsala offers training and advice to vets using ACS treatment for the first time. Please ask for more details.